

# ESI-MS Identity Confirmation Workflow

Rev. 3 · Effective 2025-01 · Quality Assurance Department

## 1. Purpose

This procedure describes the confirmation of peptide identity by electrospray ionization mass spectrometry (ESI-MS). Identity is established by agreement between the observed monoisotopic mass and the theoretical monoisotopic mass calculated from the declared sequence.

## 2. Instrument

High-resolution time-of-flight mass spectrometer operated in positive ion mode. Capillary voltage 3.5 kV, source temperature 120 °C, desolvation temperature 350 °C, desolvation gas 800 L/h.

## 3. Calibration

The instrument is calibrated daily across the  $m/z$  range 100–2000 using sodium formate clusters. Mass accuracy must be within  $\pm 2$  ppm of theoretical values for the calibration standard before sample analysis proceeds.

## 4. Sample Preparation

Dilute the HPLC working solution 1:100 in 50:50 acetonitrile:water with 0.1% formic acid to a final concentration of approximately 1  $\mu\text{g/mL}$ . Infuse at 10  $\mu\text{L/min}$ .

## 5. Acquisition

Acquire 60 seconds of TOF-MS data at 1 Hz. Combine scans, deconvolute the multiply-charged envelope, and report the monoisotopic mass of the principal species.

## 6. Acceptance Criteria

The deconvoluted monoisotopic mass must agree with the theoretical mass within  $\pm 1$  Da. For peptides above 5,000 Da, the average mass is reported and must agree within  $\pm 2$  Da. No unassigned peaks above 5% relative abundance are permitted in the deconvoluted spectrum.

## 7. Reporting

The Certificate of Analysis lists: theoretical monoisotopic mass, observed monoisotopic mass, charge states detected, and a thumbnail of the deconvoluted spectrum. Raw acquisition files are archived for seven years.

© Peptilune. This document is provided for informational purposes regarding analytical methodology and product handling.